



Gas Chromatography Analysis

Dear Customer,

Attached you will find the results of gas chromatographic analysis based on the information provided by our Producer/Laboratory for:

Sample Reference **1820**

We testify that the Chromatographic Analysis (GC) for the mentioned sample is identical with:

Product No.:	1820
Product Name:	Caraway (Carvi) organic
Botanical Name:	Carum carvi
Batch No.:	1019826

Bühl, 2020-03-25

Ayus GmbH
Am Dreschschopf 1 D-77815 Bühl/Moos
Tel.: +49-(0)7227-60099-0 Fax: +49-(0)7227-60099-99
E-Mail: info@oshadhi.eu Web: www.oshadhi.eu

Dr. David Nayan

It contains the following:

1- β - Pinene	0.106
2- D- Limonene	5.54
3- Linalool	0.10
4- Cis- Dihydrocarvone	1.11
5- Trans- Dihydrocarvone	0.37
6- Cis- Dihydrocarveol	0.40
7- Cis- Carveol	0.38
8- Trans- Carveol	0.38
9- D- Carvone	90.81