



Gas Chromatography Analysis

Dear Customer,

Attached you will find the results of gas chromatographic analysis based on the information provided by our Producer/Laboratory for:

Sample Reference **2783**

We testify that the Chromatographic Analysis (GC) for the mentioned sample is identical with:

Product No.:	2783
Product Name:	Tea-Tree, m.e. (Rosalina, linalool)
Botanical Name:	Melaleuca ericifolia
Batch No.:	1019355

Bühl, 2020-05-29


Ayus GmbH
Am Dreschschopf 1 | 77815 Bühl | Germany
Tel.: +49 7227 60099-0 | Fax: +49 7227 60099-99
E-Mail: info@ayus.team | www.ayus.team

Dr. David Nayan

α -pinene	14.26
α -terpinene	0.99
limonene	14.28
p-cymene	0.86
1,8 cineole	15.38
γ -terpinene	2.90
terpinolene	0.39
linalool	43.51
terpinen-4-ol	4.39
α -terpineol	0.34
aromadendrene	0.12
allo-aromadendrene	0.04
ledene	0.04
d-cadinene	0.08
globulol	0.02
viridiflorol	0.01